## NAViSTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter

Instructions for Use NAViSTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter	English	Page	6
Mode d'emploi Cathéter de diagnostic/d'ablation à extrémité orientable NAViSTAR® RMT DS	Français	Page	9
Gebrauchsanweisung NAViSTAR® RMT DS Diagnose-/Ablationskatheter mit lenkbarer Spitze	Deutsch	Seite	12
Istruzioni per l'uso Catetere diagnostico/ablatore a punta direzionabile NAViSTAR® RMT DS	Italiano	Pagina	15
Instrucciones de uso Catéter de diagnóstico/ablación NAViSTAR® con punta orientable RMT DS	Español	Página	18
Instruções de Utilização Cateter de Diagnóstico/Ablação com Ponta Direcionável NAViSTAR® RMT DS	Português	Página	21
Gebruiksaanwijzing NAViSTAR® RMT DS katheter met stuurbare tip voor diagnostiek en ablatie	Nederlands	Pagina	24
Brugsvejledning NAViSTAR® RMT DS diagnosticerings-/ablationskateter med styrbar spids	Dansk	Side	27
Käyttöohjeet NAViSTAR® RMT DS Ohjattavakärkinen diagnostinen-/ablaatiokatetri	Suomi	Sivu	30
Bruksanvisning NAViSTAR® RMT DS-diagnostik/ablations-kateter med rörlig spets	Svensk	Sida	33
Bruksanvisning NAViSTAR® RMT DS diagnostisk ablasjonskateter med styrbar spiss	Norsk	Side	36
Οδηγίες χρήσης Καθετήρας διάγνωσης/κατάλυσης με οδηγούμενο άκρο NAViSTAR® RMT DS	Ελληνικά	Σελίδα	39
Kullanma Talimatı NAViSTAR® RMT DS Tanısal/Ablasyon Yönlendirilebilir Uçlu Kateter	Türkçe	Sayfa	42
Инструкции по применению Диагностический/абляционный катетер с управляемым кончиком NAViSTAR® RMT DS	Русский	Стр.	45
Instrukcja obsługi Cewnik diagnostyczno-ablacyjny NAViSTAR® RMT DS o sterowalnej końcówce.	Polski	Strona	48
Pokyny k použití NAViSTAR® RMT DS pro diagnostiku/ablaci Katetr s řiditelným hrotem	Česky	Strana	51
Használati utasítás NAViSTAR® RMT DS Diagnosztikus/ablációs irányítható csúcsú katéter	Magyar	Oldal	54
Návod na použitie Diagnostický a ablačný ovládateľný hrotový katéter NAViSTAR® RMT DS	Slovensky	Strana	57
Navodila za uporabo Diagnostično/ablacijski kateter z vodljivo konico NAViSTAR® RMT DS	Slovensko	Stran	60
Указания за употреба Диагностичен/абляционен катетър с управляем връх NAViSTAR® RMT DS	Български	Страница	63
Instrucțiunile de utilizare Cateterul NAViSTAR® RMT DS cu vârf direcționabil pentru diagnosticare/ablațiune	Română	Pagina	66
Kasutusjuhend NAViSTAR® RMT DS diagnoosi/ablatsiooni juhitava otsaga kateeter	Eesti	Lk	69

Lietošanas norādījumi NAViSTAR® RMT DS diagnostikas/ekstirpācijas katetrs ar vadāmu smaili	Latviski	Lappuse	72
Naudojimo instrukcija NAViSTAR® RMT DS diagnostinis/abliācijas kateteris su valdomu galiuku	Lietuvių k.	puslapis	75
使用说明 NAViSTAR® RMT DS 诊断/消融用可调弯导管	简体中文	页码	78
사용법 NAViSTAR® RMT DS 진단/전극도자절제 조종 가능 팁 카테터	한국어	페이지	81
使用說明 NAViSTAR® RMT DS 可控式診斷/電燒導管	繁體中文	頁碼	84
Instruksi Penggunaan Kateter Diagnostik/Ablasi yang Ujungnya Dapat Dikemudikan NAViSTAR® RMT DS	Bahasa Indonesia	Halaman	87

## STERILE EO

Sterilized with ethylene oxide gas / Stérilisé à l'oxyde d'éthylène / Sterilisiert mit Ethylenoxidgas / Sterilizzato con gas ossido di etilene / Esterilizado con gas óxido de etileno / Esterilizado com gás de óxido de etileno / Gesteriliseerd met ethyleenoxidegas. / Steriliseret med ætylenoxidgas / Sterilisoitu etyleenioksidikaasulla / Steriliserad med etylenoxidgas / Steriliseret med etylenoksidgass / Αποστειρώνεται με αέριο οξείδιο του αιθυλενίου / Etilen oksit gazıyla sterilize edilmiştir / Стерилизовано этиленоксидом / Sterylizowane gazowym tlenkiem etylenu / Sterilizovano plynným etylénoxidem / Etilénoxid gázban sterilizálva / Sterilizované etylenoxidovým plynom / Sterilizirano z etilenoksidom / Стерилизиран с газ етиленов оксид / Sterilizat cu oxid de etilenă / Steriliseeritud etüleenoksiidgaasiga / Sterilizēts ar gāzveida etilēnoksidu / Sterilizuoti etileno oksido dujomis / 已采用环氧乙烷气体灭菌 / 에틸렌옥사이드 가스로 멸균 처리 / 已採用環氧乙烷氣體滅菌 / Disterilkan dengan menggunakan etilen oksida

For single use only / À usage unique / Nur zum Einmalgebrauch / Exclusivamente monouso / Para un solo uso / Destina-se apenas a uma única utilização / Uitsluitend voor eenmalig gebruik / Kun til engangsbrug. / Vain kertakäyttöön / Endast för engångsbruk / Kun til engangsbruk / Για μία χρήση μόνο / Tek kullanımlıktır / Только для одноразового использования / Instrument przeznaczony wyłącznie do jednorazowego użytku / Určeno k jednorázovému použiti / Kizárólag egyetlen használatra! / Len na jednorazové použitie / Samo za enkratno uporabo / Само за еднократна употреба / Exclusiv de unică folosință / Üksnes ühekordseks kasutamiseks / Tikai vienieizėjai lietošanai / Tik vienkartinio naudojimo / 只能使用一次 / 일회용 / 只能使用一次 / Hanya untuk digunakan satu kali

Refer to accompanying Instructions for Use / Consultez le mode d'emploi joint / Siehe zugehörige Anwendungsrichtlinien / Fare riferimento alle Istruzioni per l'uso allegate / Consulte las Instrucciones de uso adjuntas / Consultar as Instruções para Utilização incluídas / Raadpleeg de bijbehorende gebruiksaanwijzing. / Se medfølgende brugervejledning / Katso käyttöohjeita / Se medföljande bruksanvisning / Se medfølgende bruksinstruksjoner / Ανατρέξτε στις συνοδευτικές οδηγίες χρήσης / Beraberindeki Kullanma Talimatı'na bakınız / См. Инструкцию по применению, входящую в комплект поставки / Sprawdzić w załączonej Instrukcji stosowania / Viz příložený návod k použití / Használja útmutatóul a mellékelt Használati utasítást / Odkazujeme na sprievodný Návod na použitie / Upoštečajte priložena navodila za uporabo / Консультуйте се с указанията за употреба / Consultați Instrucțiunile de utilizare care însoțesc produsul / Tutvuge kaasasoleva kasutusjuhendiga / Iepazīstieties ar pievienotajiem lietošanas norādījumiem / Žr. pridedamą Naudojimo instrukciją / 请参阅随附的使用说明 / 동봉한 사용 설명서를 참조하십시오 / 請參閱隨附的使用說明。 / Lihat instruksi penggunaan yang terlampir

Do not use open or damaged packages / Ne pas utiliser si l'emballage est ouvert ou endommagé / Geöffnete oder beschädigte Verpackungen nicht verwenden / Non utilizzare confezioni aperte o danneggiate / No utilice envases abiertos o estropeados / Não utilizar embalagens abertas ou danificadas / Geopende of beschadigde verpakkingen niet gebruiken. / Bør ikke bruges, hvis emballagen har været åbnet eller er beskadiget. / Älä käytä avattuja tai vaurioituneita pakkauksia / Använd inte öppnade eller skadade förpackningar. / Bruk ikke åpne eller skadde pakninger / Μη χρησιμοποιείτε αν η συσκευασία έχει ανοίχτεί ή υποστεί βλάβη / Использование содержимого открытых или поврежденных упаковок запрещено / Açık veya hasarlı ambalajları kullanmayınız / Nie używać produktu, którego opakowanie było otwarte lub uszkodzone / Nepoužívejte, jsou-li obaly otevřené či poškozené / Ne használjon nyitott vagy sérült csomagokat / Nepoužívajte otvorené alebo poškodené balenia / Ne uporabljajte, če je ovojnina odprta ali poškodovana / Ne използвайте отворени или повредени опаковки / A nu se utiliza produsul dacă ambalajul este deschis sau deteriorat / Ärge kasutage, kui pakend on avatud või katki / Neizmantot, ja iepakojumi ir bojāti vai atvērti / Nenaudokite atidarytą arba sugadintų pakuočių / 切勿使用已打开或破损的包装 / 포장에 개방되었거나 손상된 경우 사용하지 마십시오 / 請勿使用已開啟或破損的包裝 / Jangan menggunakan kemasan yang terbuka atau yang rusak

Store in a cool, dry place / Entreposer dans un endroit sec et frais / An einem kühlen, trockenen Ort aufzubewahren / Conservare in un luogo fresco e asciutto. / Almacénese en un lugar fresco y seco / Armazenar num local fresco e seco / Bewaren op een koele, droge plaats. / Opbevares et køligt, tørt sted. / Säilytä viileässä ja kuivassa / Förvaras på en sval, torr plats / Oppbevar på et kjølig og tørt sted. / Αποθηκεύετε σε δροσερό, ξηρό μέρος / Serin ve kuru bir yerde saklayınız / Хранить в прохладном, сухом месте / Przechowywać w chłodnym, suchym miejscu / Uchovávejte na chladném, suchém místě / Száraz, hűvös helyen tárolandó. / Skladujte na chladnom a suchom mieste / Hranite v hladnem in suhem prostoru / Да се съхранява на хладно и сухо място / A se păstra la loc răcoros și uscat / Hoida jahedas ning kuivas kohas / Uzglabāt vēsā, sausā vietā / Laikykite vėsioje ir sausoje vietoje / 貯存在阴凉干燥处 / 서늘하고 건조한 곳에 보관하십시오 / 貯存在陰涼乾燥處 / Simpan di tempat yang sejuk dan kering

Use by / Utiliser avant le / Verfallsdatum / Scadenza / Fecha de caducidad / Utilizar até / Uiterste gebruiksdatum / Anvendes inden / Viimeinen käyttöpäivämäärä / Använd före / Brukes innen / Ημερομηνία λήξης / Son kullanma tarihi / Использовать до / Użyć przed / Použit do / Felhasználató: / Použitie do / Uporabno do / Срок на годност / Utilizat de: / Kölblik kuni / Izlietot līdz / Naudoti iki / 使用期限 / 유효기간 / 使用期限 / Gunakan Sebelum

## LOT

Lot No. / N° de lot / Los-Nr. / Codice lotto / Lote No. / Lote N.º / Lotnr. / Partinr. / Eränumero / Sats nr. / Lot-nr. / Ap. парт. / Lot Numarası / Номер партии / Nr partii / Č. šarže / Tételszám: / Balik č. / Serijska številka / Партиден номер / Nr. lot / Partii nr / Sērijas Nr. / Partijos Nr. / 批号 / 로트 번호 / 批號 / No. Lot

## REF

Catalog No. / N° de catalogue / Katalog-Nr. / Codice catalogo / No. de catálogo / N.º de catálogo / Catalogusnr. / Katalognr. / Kataloginúmero / Katalognr / Katalognr. / Ap. καταλόγου / Katalog Numarası / Номер по каталогу / Nr katalogowy / Katalogové číslo / Katalógusszám: / Katalógové č. / Kataloška številka / Каталоген номер / Nr. catalog / Kataloogi number / Kataloga Nr. / Katalogo Nr. / 目录号 / 카탈로그 번호 / 目錄號 / No. Katalog

Contents: 1 / Contenu : 1 / Inhalt: 1 / Indice: 1 / Contenido: 1 / Conteúdo: 1 / Inhoud: 1 / Indhold: 1 / Sisältö: 1 / Innehåll: 1 / Innhold: 1 / Περιεχόμενα: 1 / İçerik: 1 / Содержимое упаковки: 1 / Zawartość: 1 / Obsah: 1 / Tartalom: 1 / Obsah: 1 / Vsečina: 1 / Съдържание: 1 / Conținut: 1 / Sisaldab: 1 / Saturs: 1 / Sudėtis: 1 / 內含物品: 1 / 개수: 1 / 內含物品: 1 / Isi: 1

## EC REP

European Community Representative / Représentant en Europe / Vertreter in der Europäischen Gemeinschaft / Rappresentante nell'Unione Europea / Representante en la Comunidad Europea / Representante na Comunidade Europeia / Erkend vertegenwoordiger in de Europese Unie / Repræsentant i de Europæiske Fællesskaber: / Edustaja Euroopan Unionin alueella / Representant i EU / EU-representant / Αντιπρόσωπος για την Ευρωπαϊκή Κοινότητα / Avrupa Topluluğu Temsilcisi / Представитель в ЕС / Przedstawiciel we Wspólnocie Europejskiej / Zastoupení pro Evropské společenství / Európai Közösségi képviselő: / Zástupca Európskeho spoločenstva / Predstavnik v Evropski Skupnosti / Представител за Европейската общност / Reprezentant pentru Comunitatea Europeană / Volitatud esindaja Euroopa Ühenduses: / Pārstāvis Eiropas Kopienā / Atstovas Europos Bendrijoje / 歐盟代表处 / 유럽공동체 대리점 / 歐盟代表處 / Perwakilan Komunitas Eropa

Pin Connector / Connecteur à broche / Stiftverbinder / Connettore a pin / Conector de terminales / Conector de pinos / Pinconnector / Hanstik / Nastaliitin / Stiftkontakt / Hann-nålekontakt / Σύνδεσμος ακίδων / Pin Konektörü / Штекерный разъем / złącze pinowe / Kolikový konektor / Tűs csatlakozó dugó / Pinový konektor / vtični konektor / пин конектор / Conector cu pini / Pistikühendus / Kontaktu savienotājs / Šakučių jungtis / 針连接器 / 핀 커넥터 / 針連接器 / Konektor Pin

**ELECTRODES** **EA**

Electrodes / Électrodes / Elektroden / Elettrodi / Electrodos / Eléctrodos / Elektroden / Elektroder / Elektrodit / Elektroder / Elektroder / Ηλεκτρόδια / Elektrotlar / Количество электродов / Elektrody / Elektrody / Elekródák / Elekródy / Elektrode / Електроди / Electrozi / Elektroodid / Elektrodi / Elektrodoi / 电极 / 전극 / 電極 / Elektroda

**SPACING** **mm**

Spacing / Écartement / Intervall / Spaziatura / Espaciado / Espaçamento / Afstand / Afstandsstykker / Расстояние между электродами / Väli / Avstånd / Mellomrom / Απόσταση / Aralık / Odstep / Mezera / Távolság / Rozstupy / Razmik / Пространство / Distanță / Vahekaugus / Attālums / Tarpiniai / 间隔 / 간격 / 間隔 / Pengaturan Jarak

**NAV**

Navigational Catheter / Cathéter de navigation / Navigationskatheter / Catetere di navigazione / Catéter de navegación / Cateter de navegação / Navigatiekatheter / Navigerende kateter / Navigaatiokatetri / Kateter för navigering / Navigacijos kateteris / Καθετήρας πλοήγησης / Navigasyon Kateteri / Навигационный катетер / Cewnik nawigacyjny / Navigační katétr / Navigációs katéter / Navigačný katéter / Navigacijski kateter / Навигационен катетър / Cateter pentru navigare / Navigatsioonikateeter / Navigācijas katetrs / Navigacijos kateteris / 导航导管 / 내비게이션 카테터 / 導航型導管 / Kateter Navigasional

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**SINGLE USE ONLY. DO NOT RESTERILIZE.**

### DEVICE DESCRIPTION

The NAViSTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter is a steerable, multi-electrode catheter designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency current to the catheter tip electrode for ablation purposes. For mapping, the catheter is used with an RMT compatible CARTO® RMT EP Navigation System (a magnetic field location technology) and an appropriate location reference device. For ablation, the catheter is used in conjunction with a compatible radiofrequency (RF) generator and a grounding (dispersive) pad.

The NAViSTAR® RMT DS Diagnostic/Ablation Deflectable Tip catheter has an 8 mm tip with dual temperature sensors. The device has a high-torque nylon shaft with a steerable tip section containing an array of platinum electrodes, and a thermocouple temperature sensing capability. Tip deflection is controlled by the interaction of magnets placed in the catheter tip and shaft, and the magnetic field created by the external Niobe® Magnetic Navigation System (MNS). All electrodes may be used for recording and stimulation, but only the tip electrode may be used to deliver RF energy from the generator. A location sensor embedded in the tip electrode transmits location information to the RMT compatible CARTO® EP Navigation System.

The catheter interfaces with the RMT compatible CARTO® EP Navigation System and a compatible RF generator via interface cables and a Patient Interface Unit (PIU).

For use in mapping procedures, refer to the instructions for the RMT compatible CARTO® EP Navigation System.

For use in RF ablation procedures, refer to the RF generator operating instructions.

**For use in mapping and ablation procedures, refer to the instructions for the Stereotaxis Magnetic Navigation System,**

### INDICATIONS AND USAGE

The NAViSTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter, and related accessory devices are indicated for catheter-based atrial and ventricular electrophysiological mapping, and for use with a compatible RF generator in adults and children four (4) years of age and older in creating endocardial lesions during cardiac ablation procedures to treat arrhythmias.

When used with an RMT compatible CARTO® EP Navigation System, the NAViSTAR® RMT DS Diagnostic/Ablation Tip Catheter provides location information.

### CONTRAINDICATIONS

Do not use this device:

- in patients with active systemic infection;
- via the transseptal approach in patients with left atrial thrombus or myxoma, or interatrial baffle or patch;
- via the retrograde approach in patients with aortic valve replacement.

### WARNINGS AND PRECAUTIONS

- **Significant x-ray exposure can result in acute radiation injury** as well as dose-related risk for somatic and genetic effects (dose = duration of the fluoroscopic imaging X x-ray beam intensity). Take all appropriate measures to minimize x-ray exposure to both patients and clinical staff. Careful consideration should be given to the use of this device in pregnant women.
- **Ablation from within a coronary artery can cause myocardial injury and death.** Adequate fluoroscopic visualization is necessary during the transaortic approach to avoid placement of the ablation catheter in the coronary vasculature.
- **Stroke or myocardial infarction** may occur in patients undergoing **left-sided ablation procedures**. Patients should be closely monitored during the postablation period for clinical manifestations of embolic events.
- **Implantable pacemakers and implantable cardioverter/defibrillators (ICDs)** may be adversely affected by RF ablation. ICDs should be deactivated during ablation. Have temporary external sources of pacing and defibrillation available during ablation. Exercise extreme caution during ablation when in close proximity to device leads, and perform a complete analysis of the implanted device function after ablation.
- **Complete AV block** can occur when ablating septal accessory pathways or in the treatment of AVNRT. Closely monitor AV conduction during RF energy delivery and immediately terminate energy delivery if partial or complete AV block is observed. Using catheters with **distal pair electrode spacing greater than 2 mm** may increase the risk of AV nodal damage.
- **To reduce the risk of brachial plexus injury**, physicians are advised to position the patient's arms inferiorly in the normal position with the hands down by the hips as opposed to superiorly.
- Follow all warnings and precautions applicable to the Stereotaxis Magnetic Navigation System, CARTO® System and RF, Generator especially those regarding magnetic objects in the catheter laboratory. Refer to Stereotaxis Magnetic Navigation System User's Manual for user precautions.

### PRECAUTIONS

- Cardiac ablation procedures should be performed only by appropriately trained personnel in a fully equipped electrophysiology laboratory.
- Do not attempt to operate the NAViSTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter or the RF generator prior to completely reading and understanding the applicable Instructions for Use.
- The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children.
- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown.
- When radiofrequency current is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip cleaned of coagulum. When cleaning the tip electrode, be careful not to twist the tip electrode with respect to the catheter shaft; twisting may damage the tip electrode bond and loosen the tip electrode.
- Apparent low power output, high impedance reading or failure of the equipment to function correctly at normal settings may indicate faulty application of the dispersive electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication.

#### Catheter Compatibility

- The NAViSTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter is intended for use with a compatible RF generator and Biosense Webster accessories only in the Niobe® magnetic environment.
- Read and follow the dispersive electrode manufacturer's instructions for use; the use of dispersive electrodes, which meet or exceed ANSI/AAMI requirements (HF18), is recommended.

#### Handling and Sterilization

**SINGLE USE ONLY. DO NOT RE-STERILIZE.**

- This device is packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Observe "Use By" Date. Sterilized with ethylene oxide gas.

#### Disposal

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

Catheter damage may occur due to:

- autoclaving
- resterilizing
- exposure to organic solvents
- immersing proximal handle or cable connector in fluids

The sterile packaging and catheter should be inspected prior to use. If the package or the catheter appears damaged, do not use. Contact your local Biosense Webster representative.

### Environmental and EMI Precautions

- During ablation procedures, this catheter is used in conjunction with a RF generator. Electromagnetic interference (EMI) produced by the RF generator during the delivery of RF power may adversely affect the performance of other equipment. This device is a non-ionizing emitter of radiation and may cause electromagnetic interference with other devices. In order to minimize the electromagnetic interference, the generator should be placed at least 1 meter away from any other electronic device.

### Precautions During Catheter Use

- The patient should not contact grounded metal surfaces. Leakage current from any connected device to the patient must not exceed 10 microAmps ( $\mu$ A) under any circumstances. Use only isolated amplifiers, pacing equipment, and ECG equipment or patient injury or death may occur.
- Do not use excessive force to advance or withdraw the catheter. Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade.
- Always straighten the catheter tip before insertion or withdrawal of the catheter.
- Do not insert or withdraw the catheter without placing the catheter tip in neutral position.
- Use both fluoroscopy and electrograms to monitor the advancement of the catheter to the area of the endocardium under investigation to avoid vascular or cardiac damage.
- Do not immerse the electrical connector or interface cable in fluids; electrical performance could be compromised.

### Precautions During Ablation

- Do not increase power before checking for lead connection and appropriate dispersive electrode application. Effective contact between the patient and the dispersive electrode must be verified whenever the patient is repositioned.
- Do not deliver RF energy with catheter outside the target site. The RF generator can deliver significant electrical energy and may cause patient or operator injury.
- Avoid use of electrodes and probes of monitoring and stimulating devices, which could provide paths for high frequency current. Reduce the burn hazard by placing the electrodes and probes as far away as possible from the ablation site and the dispersive electrode.
- In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode cleaned of coagulum before RF current is re-applied. Use only sterile saline and gauze pad to clean the tip.
- Do not scrub or twist the tip electrode as damage may cause catheter failure or patient injury.
- Discontinue ablation immediately and replace catheter if tip temperature fails to rise during ablation.
- The temperature sensor of the catheter measures electrode tip temperature, not tissue temperature. If the generator does not display temperature, verify that the appropriate cable is plugged into the generator. If temperature still is not displayed, there may be a malfunction in the temperature sensing system, which must be corrected prior to applying RF power.

## PATIENT TREATMENT

### Antiplatelet or Anticoagulation Use

To avoid thromboemboli, intravenous heparin is used when entering the left heart during ablation, and many physicians prescribe aspirin, less often warfarin, for about 3 months afterward. No consensus yet exists about the need for short-term anticoagulation after ablation.

### Specific Patient Populations

The safety and effectiveness of cardiac ablation has not been established in:

- Asymptomatic patients;
- patients who are pregnant; or
- nursing mothers.

## PATIENT COUNSELING INFORMATION

Patients may require anticoagulation and/or antiplatelet therapy for an indefinite period based on the patient's condition.


## HOW SUPPLIED-

The NAVISTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter is provided sterile (EtO). The device is a 7 F catheter with a usable length of 130 cm, and the following features:

Tip electrode: 8 mm

Connector type: Hypertronics® interlocking

French size: 7

### Packaging

The catheter is packaged in a heat-sealed tray. The sealed tray is then placed into a pouch, forming a double sterile barrier. The sealed pouch is placed into an outer folding carton with an instruction for use (IFU) insert.

### Storage

The NAVISTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter must be stored in a cool, dry place. Storage temperature should be between 5° and 25°C (41° and 77°F).

### Shelf-Life

Product and package testing have been conducted to support the "Use By" date printed on the product labels. **DO NOT USE** after the "Use By" date.

## DIRECTIONS FOR USE

### Physician Training

Physicians must be familiar with the techniques and appropriately trained for cardiac mapping and ablation procedures. All mapping and ablation procedures must be performed in a fully equipped electrophysiology laboratory.

### Compatible RF Generators and Accessories

The NAVISTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter should be used only with a legally marketed RF generator, which has been shown to be safe and effective for cardiac ablation. Refer to the RF generator operator's manual for detailed operating instructions for RF ablation including information in choosing between temperature or power control modes.

**Specifications for Compatible RF Generators:** Refer to the RF generator manual for detailed generator operating instructions for RF catheter ablation.

**Accessories:** Use appropriate Biosense Webster accessory cables to connect the NAVISTAR® RMT DS Diagnostic/Ablation Steerable Tip Catheter to the RMT compatible CARTO® EP Navigation System and a compatible RF generator. Use an appropriate location reference device when using the catheter with the RMT compatible CARTO® EP Navigation System. Refer to the CARTO® EP Navigation System User's Manual for detailed instructions for information about catheter interface cable, location reference device and mapping procedures.

**Handling and Preparation**

- Before use, inspect the packaging. Do not use if open or damaged.
- Using aseptic technique, remove the catheter from its package and place it in a sterile working area.
- Inspect the catheter carefully for electrode integrity and overall condition. To verify compatibility between sheath and catheter, advance catheter through sheath prior to insertion.

**Operating Instructions.**

1. Create a vascular access in a large central vessel using aseptic techniques and insert the catheter.
2. Connect the catheter to the interface cables, the RMT compatible CARTO® EP Navigation System, a compatible RF generator, and standard recording equipment using the appropriate interface cables.
3. Confirm that the magnetic fields are removed.
4. Insert catheter and advance the catheter to the area of the endocardium under investigation. Use both fluoroscopy and electrograms to aid in proper positioning.
5. Use the MNS (Magnetic Navigation System) to orient the distal tip of the catheter, while advancing and retracting the catheter.
6. The catheter tip can be deflected to facilitate positioning by using the MNS to vary tip curvature. When it has been determined that the tip electrode is in stable contact with the intended ablation site, the catheter tip electrode connection must be switched from the recording equipment to the RF generator in preparation for delivery of RF current.
7. Determine that the tip electrode is in stable contact with the intended ablation site. Circuit impedance should be approximately 100 ohms upon initiation of RF current.
8. RF current may be re-applied to the same or alternate sites using the same catheter.
9. Prior to removal of the catheter, confirm that the magnetic fields are removed.

**DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY**

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

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BIOSENSE WEBSTER, the Biosense Webster logo, NAVISTAR and CARTO are trademarks of Biosense Webster, Inc.
\*The third-party trademarks used herein are trademarks of their respective owners